

# CliniScan

## Milestone 1 — Data Understanding & Preprocessing Pipeline

---

**Project Title:** CliniScan — Chest X-ray Abnormality Detection System

**Dataset Used:** VinBigData Chest X-ray Abnormalities Detection Dataset Hosted on Kaggle

### Milestone Scope

- ✓ Dataset Understanding
  - ✓ Annotation Analysis
  - ✓ DICOM Handling
  - ✓ Image Preprocessing Pipeline
  - ✓ YOLO Annotation Conversion
  - ✓ Code Implementation
- 

## PART A — DATASET UNDERSTANDING

---

### 1. Background and Motivation

Chest radiography (CXR) is one of the most widely used diagnostic imaging techniques in modern medicine. It is essential for detecting pulmonary infections, structural abnormalities, cardiovascular enlargement, pleural diseases, and oncological findings.

However, interpreting chest X-rays requires expert radiologists and is time-intensive, especially in high-volume healthcare systems. Automated detection systems can assist clinicians by:

- Highlighting suspicious regions
- Reducing diagnostic workload
- Increasing consistency
- Supporting early detection

The VinBigData dataset was designed specifically to enable research in **automated abnormality localization** rather than simple image classification.

CliniScan builds its preprocessing foundation on this dataset.

---

2. Dataset Composition

The dataset is structured into training and testing partitions.

| Component           | Description                      |
|---------------------|----------------------------------|
| Training Set        | 15,000 DICOM chest X-ray images  |
| Test Set            | 3,000 DICOM images               |
| Annotation File     | train.csv                        |
| Annotation Style    | Bounding Boxes                   |
| Radiologists        | 3 per image                      |
| Project Subset Used | 50–60 images (development phase) |

For Milestone-1, a controlled subset of 50–60 images was used to:

- Validate preprocessing pipeline
- Test annotation conversion
- Ensure computational efficiency

---

3. Project Folder Structure (As Implemented)

```
b13-cliniscan/  
|  
| └─ data/  
|   | └─ raw/  
|   |   | └─ train/ (Kaggle DICOM images)  
|   |   |  
|   |   | └─ processed/  
|   |   |   | └─ png-images/  
|   |   |   |  
|   |   |   | └─ train.csv  
|   |  
|   | └─ docs/  
|   |   | └─ Dataset_Description.pdf  
|   |   | └─ Milestone1_Report.pdf  
|  
| └─ src/
```

```
| ├── dicom_to_png.py
| ├── preprocessing.py
| └── annotation_conversion.py
|
| ├── README.md
| ├── requirements.txt
| └── LICENSE
```

This modular structure separates:

- Raw data
- Processed data
- Source code
- Documentation

This separation improves maintainability and reproducibility.

---

## 4. Understanding DICOM Format

Images are stored in **DICOM (Digital Imaging and Communications in Medicine)** format — the global standard for medical imaging.

Unlike JPEG/PNG, DICOM files contain:

- Pixel data
- Metadata
- Acquisition parameters
- Imaging modality information

### Important Technical Detail

The field `PhotometricInterpretation` can be:

- **MONOCHROME2** → Higher pixel value = Brighter
- **MONOCHROME1** → Higher pixel value = Darker (Inverted)

If inversion is not handled correctly, the lung fields appear reversed (white background, dark lungs), which can confuse a neural network.

During DICOM → PNG conversion, inversion is corrected when required.

---

## 5. Annotation File — train.csv

The train.csv file contains bounding box annotations.

Each row represents:

- One abnormality
- Annotated by one radiologist
- On one image

### Column Description

| Column     | Description                           |
|------------|---------------------------------------|
| image_id   | Unique identifier matching DICOM file |
| class_name | Abnormality name                      |
| class_id   | Numeric label (0–14)                  |
| rad_id     | Radiologist ID                        |
| x_min      | Bounding box left coordinate          |
| y_min      | Bounding box top coordinate           |
| x_max      | Bounding box right coordinate         |
| y_max      | Bounding box bottom coordinate        |

Bounding boxes are defined in **absolute pixel coordinates** relative to the original image resolution.

---

## 6. Pathology Classes

The dataset includes 14 thoracic abnormalities + 1 normal class.

| ID | Class              | Clinical Meaning           |
|----|--------------------|----------------------------|
| 0  | Aortic enlargement | Enlarged aortic silhouette |
| 1  | Atelectasis        | Lung collapse              |
| 2  | Calcification      | Calcium deposits           |

| ID | Class              | Clinical Meaning          |
|----|--------------------|---------------------------|
| 3  | Cardiomegaly       | Enlarged heart            |
| 4  | Consolidation      | Fluid-filled alveoli      |
| 5  | ILD                | Interstitial Lung Disease |
| 6  | Infiltration       | Inflammatory opacity      |
| 7  | Lung Opacity       | Increased lung density    |
| 8  | Nodule/Mass        | Suspicious lesion         |
| 9  | Other lesion       | Misc abnormality          |
| 10 | Pleural effusion   | Fluid in pleural space    |
| 11 | Pleural thickening | Thickened pleura          |
| 12 | Pneumothorax       | Collapsed lung            |
| 13 | Pulmonary fibrosis | Lung scarring             |
| 14 | No finding         | Normal scan               |

Images labeled "**No finding**" contain no bounding boxes.

---

## PART B — PREPROCESSING PIPELINE

---

### 1. Importance of Medical Image Preprocessing

Raw medical images contain variability in:

- Resolution
- Intensity scaling
- Contrast
- Scanner artifacts
- Noise patterns

Deep learning models require:

- Uniform dimensions

- Consistent pixel distribution
- Reduced noise
- Preserved anatomical detail

Therefore, a structured preprocessing pipeline was implemented.

---

## 2. Step-by-Step Pipeline

---

### Step 1 — DICOM to PNG Conversion

Tool: pydicom + OpenCV

Purpose:

- Extract pixel array
- Correct inversion (if MONOCHROME1)
- Normalize intensity
- Save as 8-bit PNG

Why convert?

- PNG is lightweight
  - Faster loading
  - Compatible with OpenCV and PyTorch
- 

### Step 2 — Grayscale Loading

Although X-rays are grayscale, OpenCV loads images in 3-channel BGR by default.

We explicitly load as:

`cv2.IMREAD_GRAYSCALE`

This:

- Reduces memory
  - Reduces computation
  - Preserves correct format
-

### **Step 3 — Resizing to 512×512**

Original resolutions vary (some > 3000×4000).

Why 512×512?

- Uniform batch training
- Computationally efficient
- Maintains sufficient anatomical detail

Interpolation used:

`cv2.INTER_AREA`

Best for downscaling.

---

### **Step 4 — Intensity Normalization (Min-Max)**

Different scanners produce different pixel ranges.

Min-Max normalization ensures:

- Consistent brightness range
  - Better convergence during training
  - Stable gradient updates
- 

### **Step 5 — CLAHE (Contrast Enhancement)**

Chest X-rays often have low local contrast.

CLAHE improves:

- Vessel visibility
- Nodule detectability
- Lung texture clarity

Parameters:

- Clip limit: 2.0
- Tile grid: 8×8

This enhances subtle abnormalities without amplifying noise excessively.

---

### Step 6 — Gaussian Blur (Denoising)

Medical imaging contains high-frequency quantum noise.

A light:

GaussianBlur (3×3)

removes noise while preserving edges.

Kernel size is small to avoid loss of lesion boundaries.

---

### Step 7 — Z-Score Normalization

After enhancement and denoising, pixel values are standardized:

$$Z = \frac{X - \mu}{\sigma}$$

Benefits:

- Centers data distribution
  - Improves neural network convergence
  - Reduces training instability
- 

### Step 8 — Annotation Conversion to YOLO Format

Original bounding boxes are in corner format.

Converted to YOLO format:

class\_id center\_x center\_y width height

Formula:

$$\begin{aligned} center_x &= (x_{min} + x_{max})/2 \\ width &= x_{max} - x_{min} \end{aligned}$$

All values normalized between 0 and 1.

Images with "No finding" → empty .txt files.

---

## 3. Why This Order Is Important



1. DICOM conversion first — raw data extraction
2. Resize before normalization — ensures consistent scale
3. Normalize before CLAHE — improves histogram quality
4. CLAHE before blur — enhance then smooth noise
5. Annotation scaling after resize — correct coordinate mapping

Changing order may distort bounding box alignment.

---

#### 4. Implementation Files

| File                     | Purpose                                    |
|--------------------------|--------------------------------------------|
| dicom_to_png.py          | Convert raw DICOM to PNG                   |
| preprocessing.py         | Apply resizing, CLAHE, blur, normalization |
| annotation_conversion.py | Convert CSV bounding boxes to YOLO         |

All scripts are modular and reusable.

---

#### 5. Output of Milestone 1

After completing this milestone:

- All selected DICOM images converted to PNG
  - Images standardized to 512×512
  - Contrast enhanced
  - Noise reduced
  - Bounding boxes converted to YOLO format
  - Clean dataset ready for object detection training
- 

#### Conclusion

Milestone-1 successfully established a structured and clinically aware preprocessing pipeline for chest X-ray abnormality detection.

By:

- Handling DICOM-specific challenges
- Standardizing image dimensions
- Enhancing local contrast
- Converting annotations properly

the dataset is now fully prepared for model development in subsequent milestones.

This forms the technical backbone of the CliniScan detection system.