

CLINISCAN PROJECT - DEEP LEARNING FOR MEDICAL IMAGING

Milestone 4

Final Evaluation, Documentation & Deployment Report

CliniScan - Chest X-ray Abnormality Detection System

Author: Adithya Jayaram

Year: 2026

Project Overview

Dataset: VinDr-CXR (15,000 images)

Framework: Streamlit (Unified Frontend & Backend)

Core Models: ResNet18 (PyTorch) + YOLOv8s (Ultralytics)

Web Application Development

- [x] Unified Web Application deployed on Streamlit Cloud
- [x] Binary Classification (Normal vs. Abnormal) via ResNet18
- [x] Localized Object Detection via YOLOv8s
- [x] Visual Explainability via Grad-CAM
- [x] Automated PDF Report Generation
- [x] Real-time inference pipeline integration

Live Application: <https://b13-cliniscan-adithya-jayaram-ppgjlzg9mpndfqyschjatz.streamlit.app>

Repository: <https://github.com/GKSJ-AI-CliniScan/B13-CliniScan/tree/Adithya-Jayaram>

CliniScan / Milestone 4 Report / Adithya Jayaram / Infosys Springboard 6.0

TABLE OF CONTENTS

PART A - Project Overview & Milestone Objectives	3
PART B - System Architecture	3
B.1 Full-Stack Architecture Overview	3
B.2 Core Processing Architecture	4
PART C - Core Inference Pipeline (Backend Logic)	5
C.1 Technology Stack	5
C.2 Pipeline Execution Flow	5
C.3 Model Integration	5
C.4 Grad-CAM Implementation	6
PART D - User Interface (Frontend Development)	6
D.1 Technology Stack	6
D.2 Core Features	7
PART E - Application Interface - Screenshots	8
PART F - Deployment	9
F.1 Streamlit Cloud Deployment	9
F.2 Environment Configuration	9
PART G - Experiments & Error Analysis	10
G.1 Hyperparameter Tuning & Augmentation	10
G.2 Error Analysis	10
PART H - Deliverables & Project Completion	11
PART I - Conclusion & Limitations	12

PART A

PROJECT OVERVIEW & MILESTONE OBJECTIVES

Milestone 4 Objective

Milestone 4 represents the final phase of the CliniScan project, transforming the trained AI models (ResNet18 and YOLOv8s) into a fully functional, publicly accessible web application. This milestone encompasses pipeline integration, explainability features, user interface design, and production deployment.

Objective	Implementation	Status
Application Server	Streamlit unified application handling both UI and inference	✓ Complete
Model Serving	ResNet18 (Classification) + YOLOv8s (Detection) loaded via PyTorch/Ultralytics	✓ Complete
Explainability	Grad-CAM implementation for ResNet18 to highlight abnormal regions	✓ Complete
User Interface	Python-based Streamlit UI with interactive elements	✓ Complete
Image Analysis	Two-stage pipeline: Classification -> Detection (if abnormal)	✓ Complete
Reporting	Automated PDF report generation containing diagnostic results	✓ Complete
Cloud Deployment	Production deployment on Streamlit Cloud	✓ Complete

PART B

SYSTEM ARCHITECTURE

B.1 Full-Stack Architecture Overview

CliniScan utilizes a streamlined, monolithic architecture powered by Streamlit. Instead of decoupled frontend and backend servers, Streamlit handles user interactions and executes backend Python logic within the same unified environment, ensuring low latency and seamless data passing.

User Browser → Streamlit Cloud Container

- **UI Layer:** Streamlit widgets (File Uploader, Buttons, Image Display)
- **Logic Layer:** State management, pre-processing, PDF generator
- **Inference Layer:** ResNet18 (Classifier) → YOLOv8s (Detector) → Grad-CAM

B.2 Core Processing Architecture

The application operates inside a containerized environment on Streamlit Cloud. It loads both AI models into memory upon startup and executes a two-stage inference pipeline based on user input.

Component	Technology	Purpose
Web Framework	Streamlit	Unified UI rendering and routing logic
Classification	PyTorch + ResNet18	Binary classification (Normal vs Abnormal)
Detection	Ultralytics YOLOv8s	Bounding box localization for abnormal cases
Interpretability	Custom Grad-CAM	Gradient-weighted activation map generation
Data Processing	OpenCV / Pillow	Image resizing, normalization, and augmentation

Deployment	Streamlit Cloud	Public hosting linked directly to GitHub repository
------------	-----------------	---

PART C

CORE INFERENCE PIPELINE (BACKEND LOGIC)

C.1 Technology Stack

Package	Role
streamlit	Core web framework and state management
torch + torchvision	Deep learning framework + ResNet18 architecture
ultralytics	YOLOv8s model loading and inference
opencv-python-headless	Image processing and bounding box annotations
Pillow	Image loading and format conversion
fpdf / reportlab	PDF report generation logic

C.2 Pipeline Execution Flow

Stage	Input	Action	Output
1. Upload	Image File	Validate and convert to required tensor shape	Processed Image Tensor
2. Classify	Image Tensor	Pass through ResNet18	"Normal" or "Abnormal"
3. Detect	Original Image	If "Abnormal", pass	Bounding Box

		to YOLOv8s	Coordinates
4. Explain	Image Tensor	Extract gradients from ResNet18	Heatmap Overlay

C.3 Model Integration

Models are loaded efficiently using Streamlit's caching mechanisms to prevent reloading weights on every user interaction.

Cached Classification model loading

@st.cache_resource

def load_classifier():

 model = resnet18(weights=None)

 model.fc = nn.Linear(model.fc.in_features, 2) # Binary classification

 model.load_state_dict(torch.load('weights/resnet18_best.pth', map_location='cpu'))

 model.eval()

 return model

Cached Detection model loading

@st.cache_resource

def load_detector():

 return YOLO('weights/yolov8s_best.pt')

C.4 Grad-CAM Implementation

Gradient-weighted Class Activation Mapping is implemented to visualize which regions of the X-ray heavily influenced the ResNet18 model's prediction.

Hook registration on ResNet18 final convolutional layer

target_layer.register_forward_hook(save_activation)

target_layer.register_backward_hook(save_gradient)

Heatmap Generation Process

1. Calculate weights from gradients

2. Apply weights to activations

3. Apply ReLU to retain only positive features

4. Normalize and apply colormap (cv2.COLORMAP_JET) overlay

PART D

USER INTERFACE (FRONTEND DEVELOPMENT)

D.1 Technology Stack

Technology	Purpose
Streamlit UI	Main layout, sidebars, interactive widgets, and layout columns
Markdown/HTML	Custom styling, text formatting, and titles within the Streamlit app
Matplotlib / Altair	Rendering data visualizations and performance metrics

D.2 Core Features

- **X-ray Upload Zone:** An interactive drag-and-drop file uploader accepting standard image formats (PNG, JPG, JPEG). Previews the uploaded image on the screen.
- **Two-Stage Analysis Trigger:** A unified action button that initiates the deep learning pipeline, managing loading states automatically.
- **Classification Results Panel:** Displays the ResNet18 binary prediction clearly (Normal vs. Abnormal) using color-coded text and confidence scores.
- **Detection Results Panel (Conditional):** If an abnormality is detected, the UI renders the X-ray with YOLOv8s bounding box annotations directly overlaid, indicating specific lesions.
- **Grad-CAM Visualization:** A side-by-side comparison showing the original X-ray alongside the Grad-CAM heatmap, providing a visual explanation for the clinician.
- **PDF Report Generator:** A custom module that compiles the patient image, classification result, detection bounding boxes, and Grad-CAM heatmap into a downloadable PDF document.

PART E

APPLICATION INTERFACE - SCREENSHOTS

(Note: Replace these placeholders with actual screenshots of your Streamlit application)

Screenshot 1 — Main Dashboard & Image Upload

(Insert image showing the Streamlit title, sidebar, and the file upload drag-and-drop zone)

Description: The main landing interface where users can upload Chest X-rays for analysis.

Screenshot 2 — Classification & Detection Results

(Insert image showing the predicted class and bounding boxes)

Description: Application displaying a positive "Abnormal" prediction alongside the YOLOv8s output highlighting the localized abnormality.

Screenshot 3 — Grad-CAM Explainability

(Insert image showing the heatmap overlay)

Description: Visual explanation panel showing the Grad-CAM heatmap, highlighting the areas the ResNet18 model focused on.

Screenshot 4 — PDF Report Download

(Insert image showing the download button and a sample PDF)

Description: The final report generation interface allowing users to download a comprehensive diagnostic summary.

PART F

DEPLOYMENT

F.1 Streamlit Cloud Deployment

The CliniScan application is deployed via Streamlit Community Cloud. This platform provides continuous integration directly from the project's GitHub repository, automatically rebuilding the application whenever new code is pushed.

Configuration	Value
Platform	Streamlit Community Cloud
Live URL	<i>(Paste your Streamlit Cloud URL here)</i>
GitHub Repo	<i>(Paste your GitHub repo URL here)</i>
Main Entry Point	app.py
Model Weights	Stored locally in repo or via external cloud storage (e.g., Google Drive / Hugging Face Hub) depending on size
Auto Deploy	Enabled – updates on every push to the main branch

F.2 Environment Configuration

A standard requirements.txt is utilized to configure the Streamlit container environment.

```
# requirements.txt snippet
streamlit==1.32.0
torch==2.1.0
torchvision==0.16.0
ultralytics==8.1.0
opencv-python-headless==4.9.0
Pillow==10.2.0
```

PART G

EXPERIMENTS & ERROR ANALYSIS

Because Streamlit operates as a single internal pipeline, testing focused on the integration of the models, hyperparameter optimization, and analyzing edge cases within the pipeline workflow.

G.1 Hyperparameter Tuning & Augmentation

Several experiments were conducted to optimize model performance before integrating them into the final deployment:

- **Learning Rate Adjustments:** Adjusted from 0.0003 to 0.0001, yielding minor performance improvements and better convergence stability.
- **Data Augmentation Strategy:** Applied flips, rotations, brightness adjustments, and CLAHE (Contrast Limited Adaptive Histogram Equalization).
 - *Outcome:* Improved recall for abnormal cases, but slightly reduced overall accuracy.
- **Alternative Architectures:** Tested DenseNet121 via transfer learning.
 - *Outcome:* Showed lower performance due to specific dataset class imbalances compared to the chosen ResNet18 architecture.

G.2 Error Analysis

Reviewing the pipeline output highlighted key challenges in automated medical imaging analysis:

Error Type	Observation	Cause / Context
False Positives	Rib shadows misclassified	Normal anatomical structures mimicking the density of pathological abnormalities.
False Negatives	Small lesions bypassed	Subtle or highly localized patterns not captured efficiently by the global classification weights.

Note: This highlights the inherent challenge of detecting subtle patterns in complex medical imaging and justifies the inclusion of the YOLOv8s detection model to force localized evaluation.

PART H

DELIVERABLES & PROJECT COMPLETION

Milestone 4 Deliverables Checklist

- [x] **Milestone 4 PDF Report** - Theory, architecture, implementation, and deployment details.
- [x] **Application Source Code** - app.py, models, and UI logic.
- [x] **Deployment Files** - requirements.txt configuration.
- [x] **Live Application URL** - (*Your Streamlit URL*)
- [x] **Classification Module** - ResNet18 integration.
- [x] **Detection Module** - YOLOv8s integration.
- [x] **Explainability Module** - Grad-CAM integration.
- [x] **PDF Report Generation** - Automated summary generation.

Complete Project Milestone Summary

Milestone	Focus Area	Key Deliverable	Status
Milestone 1	Data & Preprocessing	15,000 VinDr-CXR images prepared, Train/Test splits established	✓ Done
Milestone 2	Baseline Model Training	Initial binary classification setup and baseline bounding boxes	✓ Done
Milestone 3	Model Optimization	Trained ResNet18 and YOLOv8s with hyperparameter tuning	✓ Done
Milestone 4	App UI + Deployment	Streamlit web application and cloud hosting	✓ Done

PART I

CONCLUSION & LIMITATIONS

Conclusion

Milestone 4 successfully demonstrates the application of deep learning in medical imaging through a complete, end-to-end pipeline. The CliniScan system is now fully deployed and accessible via a web interface, transforming the trained data models into a practical decision-support tool.

The integration of ResNet18 for classification, YOLOv8s for localized detection, and Grad-CAM for visual explainability provides a comprehensive approach to automated X-ray analysis. The deployment via Streamlit Cloud highlights its practical usability, user-friendliness, and potential for further clinical development.

Final Model Performance Summary

Model	Task	Metric	Score
ResNet18	Classification	Accuracy	70%
ResNet18	Classification	ROC-AUC	0.966
ResNet18	Classification	F1-score (Abnormal)	0.73
YOLOv8s	Detection	Precision	0.5459
YOLOv8s	Detection	Recall	0.3703
YOLOv8s	Detection	mAP@50	0.4132
YOLOv8s	Detection	mAP@50-95	0.2248

Limitations & Future Work

While the classifier shows strong performance (ROC-AUC 0.96), object detection in medical contexts remains highly challenging, reflected in the moderate mAP scores. Current limitations include binary-only classification, a constrained dataset size, and a lack of formal clinical validation.

Future work will focus on expanding to multi-class classification (identifying specific diseases), improving detection accuracy with larger datasets, conducting clinical validation studies, and establishing secure integrations into existing hospital data systems.

CliniScan AI Live Application:

<https://b13-cliniscan-adithya-jayaram-ppgilzg9mpndfqyschjatz.streamlit.app>

Interface Demo:

